


EternaLift **12D HIFU**

INJECT NEW VITALITY INTO THE SKIN

BEAUTY
SHAPE
TECHNOLOGY





Brighten
skin tone



Delicate skin



Regenerated
collagen



Smooth
fine lines



Enhance
elasticity



Fat pad
reduction



Resist
gravity



Fascia
lifting

ULTRAMICRO PULSE ULTRASONIC SYSTEM

Using the Ultra-Micro Pulsed Ultrasonic System, the way instantaneous energy switching between the transverse linear acoustic pulse and the annular acoustic wave is realised. The single machine realises the combination of point, line, and surface technology mode to achieve a full mode of 360° perfect seamless full-layer anti-ageing multi-layer lifting and tightening with precise, comfortable, fast and durable effects. In MP mode, the microsecond energy transmission mode has smaller and denser focus points than other instruments, higher security, comfortable experience, shorter recovery period, more uniform and more efficient point distribution.





TT DTIRA


Multi-dimensional anti-aging effect

Like the “pyramid”

SAFE | STABLE | CONTINUOUS

MECHANICAL EFFECT+ THERMAL EFFECT+CAVITATION EFFECT

MP Technology can more accurately transmit energy to the target skin tissue layer, generating thermal effects tightening the fascia, activating the regeneration of collagen and elastic fibers, and achieving anti ageing effects such as tightening, promoting collagen regeneration, and improving skin texture without damaging surrounding tissues.



TCP IN MICRO-PULSE MODE

TCP micro-pulse mode under the micro-pulse mode realizes tiny and accurate thermal condensation point superposition, uniformly and intensively act on all areas of the target level, and the dense and scattered energy plays a more extensive role, highlighting the over all advantages.

CIRCULAR OPERATION OPTIMIZATION HAND TOOL

The annular operating hand tool integrates the world's highest new technology ultrasonic energy output mode.

- Exquisite probe, exquisite care.
- Fusion point mode & MP mode - endless output of ultrasonic energy in circular operation.

ANNULAR RUNNING PROBE

Basic nursing effect, add new operation scheme, consolidate effect and expand adaptability should be within the scope.



13.3 inch

HD smart screen


Full dynamic design

human-computer interaction

Removable plug

Full sense of technology

Patented original design



MICRO-PULSE OPERATION TECHNOLOGY&PROBE

- MP pattern (Multi-point line) [417 points/1 engine in total]
- MP pattern (Punctuate) [17 points/engine in total]
- Mechanical effect + thermal effect
- Two-way output innovation technology
- Multiple probes with different depths of care
- Multi-mode combined nursing

The advantage of multi-dimensional flat handle are targeted at the skin layer to the superficial fascia layer, fat layer, muscle layer, and a comprehensive customisation scheme based on different skin types, operating regions, and tissue levels.

CIRCULAR OPERATION OPTIMIZATION HAND TOOL&PROBE

- MP Pattern (Multi-point line) [417 points/1 engine in total]
- MP pattern (Punctuate) [12 points/1 engine in total]
- Innovative technology of micro transducer
- Endless operation of the ring
- Multiple probes with different depths of care
- Multi-mode combined nursing

The world's first annular "micro-pulse ultrasound" collaborative technology, ultrasonic linear lattice focusing, can maximise the operation effect. In addition to lifting, it also brings extra firmness; Compared with traditional ultrasound mode, multiple mode can effectively reduce pain





SPECIAL EYE PROBE

- 5.5 Mhz 2.0 mm special eye probe
- Unique ultra -thin probe technology
- World professional R&D application

Acting as an eye spot probe, the eye spot probe refines and solves problems such as shallow fine lines, deep genuine lines, dark circles, and eye bags.

FACE PROBE ACTS ON SKIN

Normal (Punctate)

Micro Plus (Linear)

- Apply gel to the skin before operation
- Operating handle vertical skin
- Do not leave the skin during operation

T- 1.5 (face) 5.5MHz


For the skin layer, superficial lines, skin color, skin texture, and pore shrinkage.

T- 3.0 (face) 5.5MHz

It directly reaches the dermis and stimulates the continuous proliferation of collagen.

T- 4.5 (face) 5.5MHz

Reach directly to the SMAS fascia layer, and tighten the muscle base to achieve tight lifting.



THE BODY PROBE IS APPLIED TO THE SKIN

Normal (Punctate)

Micro Plus (Linear)

- Apply gel to the skin before operation
- Operating handle vertical skin
- Do not leave the skin during operation

T- 6.0 (body) 2MHz

For the skin layer, skin color, skin texture, and pore shrinkage.


T- 9.0 (body) 2MHz

It directly reaches the dermis and stimulates the continuous proliferation of collagen.

T- 13 (body) 2MHz

Reach directly to the SMAS fascia layer, and tighten the muscle base to achieve tight lifting.

(5)



- Apply gel to the skin before operation
- Operating handle vertical skin
- Do not leave the skin during operation

T- 1.5 (face) 5.5MHz


For the skin layer, superficial lines, skin color, skin texture, and pore shrinkage.

T- 3.0 (face) 5.5MHz

It directly reaches the dermis and stimulates the continuous proliferation of collagen.

T- 4.5 (face) 4MHz

Reach directly to the SMAS fascia layer, and tighten the bottom of the muscle to achieve a tight lift.




PERSONAL TAILOR FULL-LAYER ANTI-AGING

According to the age group, skin status, and customer needs, point, line, and surface can be combined to achieve a full-layer, all-dimensional anti-ageing private customisation program.



- Activation of epidermis and superficial dermis
- Deep dermal activation
- Instantly tighten the fascia
- Enhance the support of fascia layer and reduce the sagging tissue
- Improve cell activity and increase cell water storage capacity
- Increase facial volume, stimulate collagen regeneration, elastic fiber mesh regeneration, plump, elastic



SINGLE MODE NURSING


MP mode (multi-point line formation)
[417 points/1 engine in total]

Dot mode can also be selected
[17 points/1 round in total]

The mode can be changed freely by touching the interface of the same probe.

MULTI-MODE COMBINED NURSING

The micro-pulse realises the superposition of tiny and precise thermal condensation points, and uniformly and intensively acts on all areas of the target level. The dense and scattered energy acts more widely, highlighting the overall advantages.



MICROPULSE OPERATION NURSING

Basic treatment efficacy, add new operation scheme, consolidate the effect and expand the scope of application.

The annular MP mode can be selected as the annular point mode, and the same probe can freely change the mode by touching the interface.



PRODUCT CONFIGURATION

Product Name : EternaLift 12D HIFU

Model : TT MPTSTL

Power Supply : AC 110V AC 220V

Power: 650W

Frequency : 2 MHz 4 MHz 5.5 MHz

Instrument weight : 10 kg

Gross product quantity : 20 kg

Package size : 20 kg

TECHNOLOGY



Instrument



Handle



Handle



Handle



Probe



Probe



Footrest



Power cord



Instructions



Packaging box

PRECAUTIONS AFTER OPERATION

- It is normal to have temporary slight redness and swelling after operation; Pay attention not to go to high temperature environments such as hot sauna, yoga and strenuous exercise within 1 week after operation, and do not expose to the sun.
- Try not to apply make-up on the day after operation (if make up needs to be removed, pay attention to the strength, and gently remove makeup oil to clean the face).
- Do not clean your face with overheated water within a week (water that does not exceed body temperature in enough); Pay attention to water supplement. It is recommended to apply water supplement facial mask every morning and evening for 3 days to 1 week; Skin care and moisturising products should be wiped 2-3 times as much as usual.
- Do not eat anything with heavy colour such as soy sauce within a week, and forbid alcohol, spicy, seafood and other strong foods.
- It is not recommended to operate similar instruments at the same place on the same day after operation.
- After operation, it is recommended to use laser and other instruments after skin repair. (The same part shall be adjusted according to the skin condition of the day).
- There may also be a small white energy rash or skin erythema at some operation sites, which are normal, mild and temporary reactions (can be treated with ice or growth factor repair products according to the specific situation).



Yolanda Cencia
Director

Ph: (07) 5689 0212
Mob: 0404 883 848
Email: yola@beautyshape.com.au


YC Empire Pty Ltd

ABN: 45 664 461 625

30109/9 Lawson St, Southport, QLD
4215, Australia



SHAPE
TECHNOLOGY



EternaLift

12D HIFU

